AUSTRIA
Bauerfeind Ges.m.b.H.
Hainburger Straße 33
1030 Wien
P +43 (0) 800 4430-130
F +43 (0) 800 4430-131
E info@bauerfeind.at

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en english

Dear patient,

Thank you for choosing a Bauerfeind product. We work to improve the medical effectiveness of our products every day – because your health is very important to us. Please read and observe these instructions for use carefully. If you have any questions, please contact your physician or medical retailer.

Intended purpose
The GenuLoc is a medical device. It is an orthosis[1] for immobilizing the knee joint – either with the leg extended or at a preset flexion angle.

⚠ **Fitting of the GenuLoc may only be carried out by a professional[2], otherwise product liability cannot be assumed.**

Indications
- Pre-operative, post-operative, post-traumatic immobilization
- Non-surgical treatment of patella fractures

Notes before first usage
If you wish to use your Bauerfeind product for acute symptoms, injuries or healed previous conditions, you urgently need to obtain professional medical advice and observe it before using the product. Restrictions may apply to you which you must take into account and which the professional[1] will clarify.
In case of identified or suspected problems with blood pressure and intense physical exercise, prolonged use of the product can have additional effects on your blood pressure.
After you have been prescribed/you have purchased the product, only use it according to the indication and under observance of any additional instructions given by your physician/professional[1]. If you wish to use the product in combination with medication, other medical devices or other treatments, consult your physician/professional[1] beforehand.
If any additional conditions or injuries occur during the period of application of the product, also consult your physician immediately.

Risks of using this product
⚠ **Important information**
- Correct usage / donning of the product is assumed. Any therapeutic aids – supports and orthoses – applied externally to the body can, if tightened excessively, lead to local pressure symptoms or, rarely, constrict blood vessels or nerves. If you experience this, it can be helpful to loosen the straps slightly and immediately consult your physician and professional[2].
- The knee orthosis GenuLoc is a prescribed product that should be applied under a physician's guidance. GenuLoc should only be worn in accordance with these instructions for use and for the listed areas of application.
- Product liability cannot be assumed in the event of inappropriate use.
- No improper modifications may be made to the product. Failure to comply with this requirement may adversely affect product performance, thereby excluding any product liability.
- If you notice any changes or an increase in symptoms while wearing the product, stop any further use and contact your physician.
 Caution: When loosening the straps or tensioning and fastening systems, or when doffing the product or adjusting configured limitations of movement, there is not sufficient support / sufficient protection of the body part requiring treatment. In this case, particularly avoid placing any strain on this body part.
- If you would like to use your Bauerfeind product due to acute symptoms / injuries, make sure to obtain and follow the advice of a professional before first using the product. Restrictions may exist given the specific indication which the professional can clarify. In particular, the ability to operate vehicles, other means of transport or machinery may be limited. In case of doubt, we recommend temporarily refraining from the activities mentioned above.
- Combining the support with other products (e.g. compression stockings) must be discussed beforehand with your physician.
- Do not allow the product to come into contact with ointments, lotions, or substances containing grease or acids.

Contraindications
Side effects involving the body as a whole have not been reported to date.
- No clinically significant hypersensitivity reactions have been reported to date. Generally, treatment with a knee immobilization splint is prescribed by a physician. If this has not occurred and you have any of the following conditions, consult your physician before donning and wearing the product:
 – Impaired sensation and circulatory disorders in the legs, e.g. diabetes mellitus
 – Skin disorders / injuries in the relevant part of the body, particularly if inflammation is present. In case of any raised scars with swelling, redness or excessive heat build-up
 – Severe varicose veins (varicosis), particularly with impaired return flow of blood to the heart
 – Impaired lymphatic drainage – including unspecified soft tissue swelling located in distal direction of the fitted aid

Application instructions
GenuLoc is available in one universal size and can be adjusted individually. For this reason, size determination is not carried out.

Before donning the orthosis for the first time:
(only by a professional[2])
- Before using the product, the professional[2] must check for the need of thrombosis prophylaxis.
- The preset bending angle must be checked by the professional[2] and adjusted if necessary. Please note that only the two aluminum splints in the hollow of the knee can be adjusted.
- Before donning the system, the circumference of the orthosis is adjusted using the lateral Velcro straps. The side parts must be attached in the same way on both sides (Image 1 / 2).
- If the use of a pad behind the knee is required (please ask your physician about this), attach this to the Velcro area provided on the inside of the orthosis before donning the product.

Donning of the orthosis
- When positioning the system on the leg, make sure that the distance from the middle straps to the middle of the knee caps is uniform (Image 3 / 4).
- To close the system, start with the lower strap near the joint (Image 5).
- Then close the upper strap near the joint (Image 6).
- To don the orthosis completely, close the lower strap away from the joint and then the upper strap away from the joint (Image 7 / 8).
- To firmly fix the orthosis in place, first evenly pull tight the two lower straps running in opposite directions. Proceed the same way with the two upper straps (Image 9).

Doffing of the orthosis
- When doffing the GenuLoc, the straps are opened in the reverse order (Image 9, 8, 7, 6, 5). For your own safety, carry this out while sitting.

Cleaning instructions
- Please only wash the orthosis by hand at 30 °C using a mild detergent. All components can stay on the orthosis during washing. Leave the orthosis to air dry.
- Regular care guarantees the ideal effectiveness of the product. Also observe the instructions on the label behind the pad at the back of the knee.
- Never expose the product to direct heat / cold as this may impair its effectiveness.
- If you have any complaints, please contact your distributor. Only clean goods are handled.

Part of the body this product is used for
In accordance with the indications (knee). Please see intended purpose.

Note on maintenance
If handled and cared for correctly, the product is practically maintenance-free.

Technical specifications / parameters
GenuLoc is a stabilizing orthosis for immobilization of the knee joint. It consists of supporting material, four stabilizing stays and four straps.

Notes on reuse
This product is intended to be fitted to one patient only.

Warranty
The legal provisions of the country in which the product was purchased apply. If you suspect a warranty claim, please first contact the person from whom you purchased the product directly. The product must be cleaned before making a warranty claim. If the instructions for handling and caring for the product have not been observed, the warranty may also be affected or excluded.
The warranty does not cover:
- Off-label use
- Failure to observe the instructions of the physician/professional[1]
- Arbitrary product modifications

Disclaimer
In general, do not self-diagnose, self-treat or self-medicate unless you are a professional[1].
Before first using our medical device, actively consult a physician or professional[1] since this is the only way to assess the individual effect of the product on your body and determine any possible risks of application/side effects. Please observe the advice of this professional[1] as well as all notes and instructions in this document or its online version (including texts, images, graphics, etc.), including any excerpts.

Obligation to report
Due to regional legal regulations, you are obliged to report any serious incidents involving the use of this medical device to both the manufacturer and the competent authority for medical devices (in Germany: Bundesinstitut für Arzneimittel und Medizinprodukte - Federal Institute for Drugs and Medical Devices).

Disposal
Once you have finished using the product, please dispose of it according to local specifications.

Material content
Polyamide (PA), Aluminum (AlMgSi1), Polyurethane (PUR), Polyoxymethylene (POM), Cotton (CO), Polyester (PES), Elastane (EL)

MD – Medical device
UDI – DataMatrix Unique Device Identifier UDI

Date of information: 2026-09

[1] Orthosis = orthopedic aid used to stabilize, relieve, immobilize, control, or correct the limbs or the torso.

[2] A professional is any person who is authorized according to the state or institutional regulations as well as other provisions relevant to them for fitting and instruction in the use of orthopedic aids. As the manufacturer, we strongly recommend that this person be trained in orthopedics or possess comparable acquired competencies as well as dispose of orthopedic infrastructure.

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