MEDICAL
BRAND
SINCE 1929

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Stand der Information: 2026-07

Dear patient,
Thank you for choosing a Bauerfeind product.
We wish you a speedy recovery.

Please read and observe these instructions for use closely and in their entirety. If there is something you do not understand or if you have any questions, do not use the product until you have consulted your physician or professional[1] again. Keep these instructions for use in case you want to look something up later.

Intended purpose
The GenuTrain P3 is a medical device. It is a support for centering and securing the patella with an adjustable corrective strap.

Indications
- Knee joint irritation (patellofemoral pain syndrome, patellar chondropathy, chondromalacia patellae, patella tip syndrome/jumper's knee, runner's knee)
- Retropatellar osteoarthritis
- Patellar lateralization
- Patellar subluxation/patellar dislocation (post-traumatic)
- High-riding patella (patella alta)
- Post-operative care (e.g. following stabilization of the patella)
- Feeling of instability

Contraindications
No clinically significant hypersensitive reactions have been reported to date. This aid should be donned and worn only after consulting your physician if you have one of the following medical conditions:
- Skin disorders/injuries in the relevant part of the body, particularly if inflammation is present. In case of any raised scars with swelling, redness or excessive heat build-up
- Varicose veins
- Impaired sensation in the leg due to neurological abnormalities
- Impaired circulation in the leg (e.g. diabetes mellitus)
- Impaired lymphatic drainage and soft tissue swelling of unknown origi

Notes before first usage
If you wish to use your Bauerfeind product for acute symptoms, injuries or healed previous conditions, you urgently need to obtain professional medical advice and observe it before using the product. Restrictions may apply to you which you must take into account and which the professional[1] will clarify.
In case of identified or suspected problems with blood pressure and intense physical exercise, prolonged use of the product can have additional effects on your blood pressure.
After you have been prescribed/you have purchased the product, only use it according to the indication and under observance of any additional instructions given by your physician/professional[1]. If you wish to use the product in combination with medication, other medical devices or other treatments, consult your physician/professional[1] beforehand.
If any additional conditions or injuries occur during the period of application of the product, also consult your physician immediately.

Risks of application ⚠
- Do not arbitrarily modify the product as this may cause it to not help as expected and/or worsen your health.
- The product must not be applied too firmly as this may cause local pressure symptoms or, in rare cases, constrict underlying blood vessels or nerves.
- If you notice any changes or an increase in symptoms related to the application of the product (such as tissue congestion, skin irritation, swelling, etc.), pause the application and consult your physician/professional[1].
- If you completely doff the product, treatment is insufficient and additional damage to the part of the body affected may occur during additional strain. Be especially careful when this is the case.
- Additionally, restrictions due to the indication may apply which the professional[1] will explain to you. In particular, driving vehicles, other means of transportation or machines may only be possible under restrictions. The same applies to the use of ladders, tools and sports equipment. When in doubt, we recommend to suspend these activities.
- Persons requiring supervision due to their age or mental/physical impairment must be supervised during every application of the product.
- Avoid contact of the product with fatty and acidic agents, ointments, lotions, bases, sand or salt water.
- If hypersensitivity has been observed with you, consult your physician before the first application of the product. If hypersensitivity/an allergic reaction occurs in relation with the application of the product, pause the application and consult a physician immediately.
- Never expose your product to direct sunlight, heat or cold for longer periods of time.

Application instructions
- Observe the instructions for donning and doffing of the product as well as any additional instructions for use given by your physician/professional[1] under any circumstances.
- Since the product unfolds its effects especially during physical activity, doff it during longer rest phases.
- Check the product for integrity before every use. Only intact products may be used.
- Make sure to have your physician/professional[1] show you how to don/doff the product and explain to you how it works until no doubts remain. If the affected part of the body changes in diameter, the product may have to be readjusted/newly prescribed.

Donning of the product
① Grip the support by the donning aids at the top ends of the lateral stays on both sides ② and pull the support up over the knee. ③ Position the ring-shaped pad in such a way that it encloses the kneecap without exerting pressure. ④ If required, correct the fit of the support.

Securing the kneecap in the desired position
⑤ To ensure that the kneecap is in the perfect position, bend the leg slightly (max. 20°). Release the tensioning strap of the corrective strap from the Velcro area on the inside of the support. Pull the strap diagonally upward until your kneecap is in the desired position that is comfortable for you. ⑥ Refasten the end of the strap onto the Velcro area.

Doffing of the product
⑦ Grip the support by the aid for taking off your stockings at the bottom ends of the stays on both sides and pull the support downwards. Make sure not to turn the support inside-out.

Cleaning instructions
Observe the supplementary instructions on the sewn-in label/sticker on your product. Wash the product at 30 °C/86 °F. Regular care ensures the ideal effectiveness of the product.

Area of use
Knee joint

Maintenance instructions
If handled and cared for correctly, the product is practically maintenance-free.

Assembly and fitting instructions
The product is supplied ready for use

Technical specifications/parameters
The GenuTrain P3 consists of a knitted section with an integrated massage pad around the kneecap, a corrective strap on top of this and an oval massage pad at the upper edge of the support. Spare parts are not available.

Notes on reuse
The product is intended for the individual treatment of one patient only.

Warranty
The legal provisions of the country in which the product was purchased apply. If you suspect a warranty claim, please first contact the person from whom you purchased the product directly. The product must be cleaned before making a warranty claim. If the instructions for handling and caring for the product have not been observed, the warranty may also be affected or excluded.

The warranty does not cover:
- Off-label use
- Failure to observe the instructions of the physician/professional[1]
- Arbitrary product modifications

Disclaimer
In general, do not self-diagnose, self-treat or self-medicate unless you are a professional[1].
Before first using our medical device, actively consult a physician or professional[1] since this is the only way to assess the individual effect of the product on your body and determine any possible risks of application/side effects. Please observe the advice of this professional[1] as well as all notes and instructions in this document or its online version (including texts, images, graphics, etc.), including any excerpts.

Obligation to report
Due to regional legal regulations, you are obliged to report any serious incidents involving the use of this medical device to both the manufacturer and the competent authority for medical devices (in Germany: Bundesinstitut für Arzneimittel und Medizinprodukte - Federal Institute for Drugs and Medical Devices).

Disposal
Please dispose of the product in accordance with local regulations.

Material content
Polyamide (PA), Silicone (SI), Thermoplastic polyurethane (TPU), Elastane (EL), Thermoplastic elastomer (TPE-SEPS), Polyurethane (PUR), Cotton (CO), Polydimethylsiloxane-polyurea block copolymer (based on PDMS-urea copolymer), Polyester (PES)

MD – Medical Device
UDI – DataMatrix identifier as UDI

1 A professional is any person who is authorized to adjust supports and orthoses and instruct others in their use in accordance with the applicable state and institutional regulations and other guidelines. As the manufacturer, we strongly recommend that this person disposes of orthopedic technical training or comparable acquired skills, as well as orthopedic technical infrastructure.

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