ENGLISH

BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

INTENDED USER
The intended user should be a licensed healthcare professional, the patient or the patient's caregiver. The user should be able to read, understand and be physically capable of performing the directions, warnings and cautions in the information for use.

INTENDED USE/INDICATIONS
The Donjoy Bellystrap is designed to provide comfort and support to the abdomen. It may be suitable for strains or muscle weakness of the abdominal area and for post-operative abdominal surgery.
- Strains or muscle weakness of the abdominal area
- Post-operative abdominal surgery

PERFORMANCE CHARACTERISTICS
- Providing support and stabilization to the limb or body segment.
- Providing partial immobilization or controlled movement to the limb or body segment.

CONTRAINDICATIONS
- If you have skin disorders/injuries in the relevant part of the body, the product should be fitted and worn only after consultation with a medical professional.
- Do not use if you are allergic to any of the materials contained within this product.

WARNING AND PRECAUTIONS
- We recommend that the first application be made under supervision of a healthcare professional.
- If you have skin conditions/injuries in the relevant part of the body, the product should be fitted and worn only after consultation with a healthcare professional.
- If you experience any discomfort, pain, swelling or sensation changes, or any unusual reactions while using this product, consult your healthcare professional immediately.
- This product is for single patient multiple use only.
- Do not use this device if it was damaged and/or packaging has been opened.
- The product should be removed during long periods of rest.
- Do not use during showering or bathing.
- Keep the product away from children.
- Do not use on open wound(s) and/or sensitive skin.
- Do not use this product if it is damaged or defective.
- Do not modify the device.
- Do not use this product in any way that may damage its integrity or use the device other than intended.
- The frequency and duration of use should be determined by a healthcare professional.

NOTE: Contact manufacturer and competent authority in case of a serious incident arising due to usage of this device.

APPLICATION INSTRUCTIONS
If possible, the abdominal support should be put on when standing.
Open the belt.
1. Place the belt around the waist. Make sure that the abdominal belt is correctly centred over your spine (the back stays should be positioned on either side of the spine). Position the front panel with rounded shape over your abdomen.
2. While maintaining the front panel in position, pull the part with the hook-and-loop fasteners and attach it to the abdominal support.
3. Adjust the hook-and-loop closure tab until desired support is obtained.
4. All the stays can be removed through the slots.
5. The level of compression can be readjusted as needed during the day.

DISPOSAL
Dispose of the device and its packaging in accordance with local or national environmental regulations.

CLEANING INSTRUCTIONS
Machine washable at gentle cycle (max 40°C) with a mild detergent.
Hook-and-loop fasteners of the straps should be closed for washing to avoid any damage.
Do not bleach. Rinse thoroughly. Air dry. Do not machine dry. Do not dry clean. Do not iron.

COMPOSITION
Polyester, Cotton, Polyamide, Steel, Polethylene, Elastane, Polyurethane

NOT MADE WITH NATURAL RUBBER LATEX

FOR SINGLE PATIENT USE ONLY.

WARRANTY DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale. To the extent the terms of this warranty are inconsistent with local regulations, the provisions of such local regulations will apply.

NOTICE: WHILE EVERY EFFORT HAS BEEN MADE IN STATE-OF-THE-ART TECHNIQUES TO OBTAIN THE MAXIMUM COMPATIBILITY OF FUNCTION, STRENGTH, DURABILITY AND COMFORT, THERE IS NO GUARANTEE THAT INJURY WILL BE PREVENTED THROUGH THE USE OF THIS PRODUCT.

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©2024 DJO, LLC