ENGLISH

Before using the device, please read the following instructions completely and carefully. Correct application is vital to the proper functioning of the device.

INTENDED USER PROFILE

- The intended user should be a licensed medical professional, the patient, the patient's caretaker, or a family member providing assistance.
- The user should be able to read, understand and be physically capable of performing all the directions, warnings and cautions provided in the information for use.

INTENDED USE/INDICATIONS

- Support and compression following acute and chronic knee injury and mild knee osteoarthritis
- Proprioceptive and neuromuscular stimulation through compression (to support balance conditions and movements)
- Prevention of injury

Provide strong compression (20-36 mmHg, French Standard ASQUAL) of limb or body segment through elastic construction.
Soft-good/semi-rigid designed to restrict motion through elastic or semi-rigid construction.

CONTRAINDICATIONS

- If you have skin disorders/injuries in the relevant part of the body, the product should be fitted and worn only after consultation with a medical professional.
- Do not use if you are allergic to materials in this product or materials containing Neoprene especially phlebitis patients.
- Impaired blood flow in the leg (e.g. with diabetes mellitus).
- Impaired lymphatic drainage or soft tissue swelling of unknown origin.

WARNING & PRECAUTIONS

- The product is intended for single patient use only.
- The product is designed to be worn directly on the skin.
- The product should not be used for extended periods of time, and removed during long periods of rest.
- Do not apply any lotions or ointments on the relevant part of the body while wearing the support.
- If pain, swelling, changes in sensation or other unusual reactions occur while using this product, you should contact a medical professional immediately.
- Do not use over open wounds.
- Do not use this device if it was damaged and/or packaging has been opened.

Note: Contact manufacturer and competent authority in case of a serious incident arising due to usage of this device.

APPLICATION INSTRUCTIONS

1. Grip edges of the support at the larger opening.
2. Pull the support over the knee as shown on the picture. (Fig. 1).
3. Adjust the fit of the support as necessary.

CLEANING INSTRUCTIONS

Machine washable (max 40°C) with a mild detergent. Do not bleach. Rinse thoroughly. Air dry. Do not machine dry. Do not dry clean. Do not iron.

COMPOSITION

Elastane, Polyamide, Polyester, Polyolefin

NOT MADE WITH NATURAL RUBBER LATEX

WARRANTY DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale. To the extent the terms of this warranty are inconsistent with local regulations, the provisions of such local regulations will apply.

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©2006 DJO, LLC
DJO LLC
DJO FRANCE SAS