DONJOY®

ENGLISH

BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

INTENDED USER PROFILE:
The intended user should be a licensed medical professional, the patient, the patient's caretaker, or a family member providing assistance. The user should be able to: Read, understand and be physically capable to perform all the directions, warnings and cautions provided in the information for use.

INTENDED USE/INDICATIONS:
The DonJoy Immo 3V brace is designed to provide support and immobilization to the knee. It may be suitable for use following a knee injury or soft tissue injury. Providing immobilization or controlled movement of the limb or body segment.

CONTRAINDICATIONS:
Do not use if you are allergic to any of the materials contained in this product.

WARNINGS AND PRECAUTIONS:
- It is recommended to have the device fitted by a suitably qualified healthcare professional.
- If pain, swelling, changes in sensation or other unusual reactions occur while using this product, you should contact your doctor immediately.
- Do not use over open wounds.
- Do not use this device if it was damaged and/or packaging has been opened.

NOTE: Contact manufacturer and competent authority in case of a serious incident arising due to usage of this device.

APPLICATION INSTRUCTIONS:
1) Open up the closure straps.

2) Position the 3 panels so as to best cover the patient's leg, placing the patellar opening over the patella. Try the immobilizer over the patient's leg.

3) Close the immobilizer and tighten the closure straps. Tighten the straps according to the circumference of the thigh and calf, making sure they are neither too tight nor too loose.

Note: Shaping the posterior stay helps improve the stability of the immobilizer.

CARE:
Hand wash in cold water with mild soap, rinse thoroughly. Air dry. Do not machine dry or bleach.

COMPOSITION:
Polyester, Polyamide, Polyurethane, Polypropylene, Aluminium.

FOR SINGLE PATIENT USE ONLY.
NOT MADE WITH NATURAL RUBBER LATEX.

WARRANTY DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale.

NOTICE: WHILE EVERY EFFORT HAS BEEN MADE IN STATE-OF-THE-ART TECHNIQUES TO OBTAIN THE MAXIMUM COMPATIBILITY OF FUNCTION, STRENGTH, DURABILITY AND COMFORT, THERE IS NO GUARANTEE THAT INJURY WILL BE PREVENTED THROUGH THE USE OF THIS PRODUCT.

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©2021 DJO, LLC
DJO FRANCE S.A.S