ENGLISH

BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

INTENDED USER PROFILE
The intended user should be the patient, a licensed medical professional, or the patient's caregiver.
The user should be able to read, understand and be physically capable of performing the directions, warnings and cautions in the information for use.

INTENDED USE/INDICATIONS
- Support and compression following injury to the patellar tendon
- Provides proprioceptive and neuromuscular stimulation through compression (to support balance conditions and movements)
- Helps in the management of overuse injuries such as Osgood-Schlatter disease and patellofemoral pain syndrome
- Patellar tendinitis
- Patellofemoral pain syndrome
- Patellar maltracking

Soft-good designed to provide mild protection and strong compression (20-36 mmHg, French Standard ASQUAL) through elastic construction.

CONTRAINDICATIONS
- Any known allergies to the medical device materials.
- Peripheral arterial diseases.
- Impaired blood flow in the leg (e.g. with diabetes mellitus).
- Impaired lymphatic drainage or soft tissue swelling of unknown origin.

WARNINGS AND PRECAUTIONS
- If you have skin conditions/trauma in the relevant part of the body, the product should be fitted and worn only after consultation with a medical professional.
- Do not apply on open wounds.
- The product is designed to be worn directly on the skin.
- Do not apply any lotions or ointments on the relevant part of the body while wearing the support.
- Remove the support during longer periods of rest.
- If pain, swelling, changes in sensation or other unusual reactions occur while using this product, you should contact your doctor immediately.
- Handle with care: avoid fastening the strap's hook closure onto the knitted sleeve as it may cause damage to the appearance of the product when strap is removed.
- Do not use this device if it is damaged and/or packaging has been opened.

NOTE: Contact manufacturer and competent authority in case of a serious incident arising due to usage of this device.

APPLICATION INSTRUCTIONS
1) Unfasten and loosen the strap.
2) Grip edges of the support at the larger opening and pull the support over the knee.
3) Adjust the fit of the support as necessary.
Ensure that the big pad encloses the kneecap and the small pad (tendon pad) is located below the kneecap over the patellar tendon.
4) Tighten the strap to set the pressure and then close it.
The strap should run over the small pad.

CLEANING INSTRUCTIONS
Machine washable at gentle cycle (max 40°C) with a mild detergent.
Hook and loop should be closed for washing.
Do not bleach. Rinse thoroughly. Air dry. Do not machine dry. Do not dry clean. Do not iron.

COMPOSITION
Polyamide, Elastane, Polyester, Steel, Polyethylene, Elastomer, Polyurethane, Cotton, Polyolefin, Silicone

NOT MADE WITH NATURAL RUBBER LATEX

WARRANTY DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale. To the extent the terms of this warranty are inconsistent with local regulations, the provisions of such local regulations will apply.

NOTICE: WHILE EVERY EFFORT HAS BEEN MADE IN STATE-OF-THE-ART TECHNIQUES TO OBTAIN THE MAXIMUM COMPATIBILITY OF FUNCTION, STRENGTH, DURABILITY AND COMFORT, THERE IS NO GUARANTEE THAT INJURY WILL BE PREVENTED THROUGH THE USE OF THIS PRODUCT.

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©2022 DJO, LLC 13-4431 REV A - 2022-01-26
DJO France SAS
MDSS GmbH