ENGLISH

BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

INTENDED USER
The intended user should be a licensed healthcare professional, the patient or the patient's caregiver. The user should be able to read, understand and be physically capable of performing the directions, warnings and cautions in the information for use.

INTENDED USE/INDICATIONS
The Donjoy Sikinstrap is designed to provide comfort and support of the lumbar spine. It may be suitable for the symptomatic relief of acute and chronic low back pain. Provides compression without significant restriction of mobility.
- Acute lower back pain
- Chronic low back pain
- Muscular weakness

PERFORMANCE CHARACTERISTICS
- Providing support and stabilization to the limb or body segment.
- Providing partial immobilization or controlled movement to the limb or body segment.
- Providing mild protection of a limb or body segment for acute and prophylactic care.

CONTRAINDICATIONS
- If you have skin disorders/injuries in the relevant part of the body, the product should be fitted and worn only after consultation with a medical professional.
- Do not use if you are allergic to any of the materials contained within this product.

WARNING AND PRECAUTIONS
- We recommend that the first application be made under supervision of a healthcare professional.
- If you have skin conditions/injuries in the relevant part of the body, the product should be fitted and worn only after consultation with a healthcare professional.
- If you experience any discomfort, pain, swelling or sensation changes, or any unusual reactions while using this product, consult your healthcare professional immediately.
- This product is for single patient multiple use only.
- Do not use this device if it was damaged and/or packaging has been opened.
- The product should be removed during long periods of rest.
- Do not use during showering or bathing.
- Keep the product away from children.
- Do not use on open wound(s) and/or sensitive skin.
- Do not use this product if it is damaged or defective.
- Do not modify the device.
- Do not use this product in any way that may damage its integrity or use the device other than intended.
- The frequency and duration of use should be determined by a healthcare professional.

NOTE: Contact manufacturer and competent authority in case of a serious incident arising due to usage of this device.

APPLICATION INSTRUCTIONS
If possible, the back support should be put on when standing.
Open the belt.
1. Place the belt around the waist, placing your hands in the finger loops located on the front closures. Make sure that the lumbar belt is correctly centred over your spine (the stays should be positioned on either side of the spine). Position the bottom part of the belt at the top of the gluteal cleft.
2. Pull the finger loops to extend the main fabric and wrap it around your waist.
3a/b. Fasten the belt with the right side over the left.
4. The 2 front stays can be removed through the slots.
5. The level of compression can be readjusted as needed during the day.

ACCESSORIES
The Donjoy Sikstrap has been designed to be used with the following optional accessories:
- Massaging back pad: The dots should face the back's curvature. Fasten the pad to the back of the belt using the hook-and-loop fasteners (arrow symbol pointed upwards).
- Shoulder straps: These allow you to open the belt when support is not needed without having to remove the entire belt. Attach the shoulder straps to the belt using the hook-and-loop fasteners.
- Hot/cold pack with pouch. Fasten the pouch to the back belt using the hook-and-loop fasteners.

DISPOSAL
Dispose of the device and its packaging in accordance with local or national environmental regulations.

CLEANING INSTRUCTIONS
Machine washable at gentle cycle (max 40°C) with a mild detergent.
Hook-and-loop fasteners of the straps should be closed for washing to avoid any damage.
Do not bleach. Rinse thoroughly. Air dry. Do not machine dry. Do not dry clean. Do not iron.

COMPOSITION
Front Closure: Polyamid, Polyester, Polypropylene, Polyurethane
Belt: Polyamide, Polyester, Polyethelene, Steel

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© 2024 DJO, LLC
DJO, LLC