ENGLISH

Before using the device, please read the following instructions completely and carefully. Correct application is vital to the proper functioning of the device.

INTENDED USER PROFILE

The intended user should be a licensed medical professional, the patient, the patient's caretaker, or a family member providing assistance.
The user should be able to read, understand and be physically capable of performing all the directions, warnings and cautions provided in the information for use.

INTENDED USE/INDICATIONS

- Support and compression following acute and chronic knee injury and mild knee osteoarthritis
- Proprioceptive and neuromuscular stimulation through compression (to support balance conditions and movements)
- Prevention of injury
- Mild medial (inner) and lateral (outer) ligament sprains of the knee
- Mild knee instability

Provide strong compression (20-36 mmHg, French Standard ASQUAL) of limb or body segment through elastic construction.
Soft-good/semi-rigid designed to restrict motion through elastic or semi-rigid construction.
Providing mild protection of a limb or body segment for acute and prophylactic care.

CONTRAINDICATIONS

- If you have skin disorders/injuries in the relevant part of the body, the product should be fitted and worn only after consultation with a medical professional.
- Any known allergies to the medical device materials.
- Peripheral arterial diseases.
- Impaired blood flow in the leg (e.g. with diabetes mellitus).
- Impaired lymphatic drainage or soft tissue swelling of unknown origin.

WARNING & PRECAUTIONS

- The product is intended for single patient use only.
- We recommend that the first application be made under supervision of a health care professional.
- The product is designed to be worn directly on the skin.
- The product should not be used for extended periods of time, and removed during long periods of rest.
- Do not apply any lotions or ointments on the relevant part of the body while wearing the support.
- If pain, swelling, changes in sensation or other unusual reactions occur while using this product, you should contact a medical professional immediately.
- Do not use over open wounds.
- Do not use this device if it was damaged and/or packaging has been opened.
- Handle with care; avoid fastening the strap's hook closure onto the knitted sleeve as it may cause damage to the appearance of the product when strap is removed.

Note: Contact manufacturer and competent authority in case of a serious incident arising due to usage of this device.

APPLICATION INSTRUCTIONS

1. Grip edges of the support at the larger opening and pull the knee support over the knee then adjust the fit of the support, as necessary.
Make sure that the circular pad encloses the kneecap.
2. Fix the circle hook of the strap to the matching velvet circle on the knee support.
When fixing the circle hook, place the strap on the opposite side of the injured ligament.
(i.e. toward the outer side of the knee for injuries of the medial (inner) ligament or toward the inner side of the knee for injuries of the lateral (outer) ligament).
3. Wrap the strap around the leg (below the knee).
4. Pull the strap behind the knee and the thigh.
5. Wrap it around the thigh.
6A & 6B. Pull the strap over the injured ligament and wrap it around the leg (below the knee).
7. Pull the strap again over the injured ligament.
Attach the hook and loop tab (with DonJoy logo) to the strap itself on the thigh.

CLEANING INSTRUCTIONS

Machine washable (max 40°C) with a mild detergent.
Hook and loop should be closed for washing.
Do not bleach. Rinse thoroughly. Air dry. Do not machine dry. Do not dry clean. Do not iron.

COMPOSITION

Elastane, Polyamide, Cotton, Polyester, Polyolefin, SBR Foam, Polyurethane, Thermoplastic Polyurethane, Elastomer, Polypropylene, Steel, Silicone

NOT MADE WITH NATURAL RUBBER LATEX

WARRANTY DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale. To the extent the terms of this warranty are inconsistent with local regulations, the provisions of such local regulations will apply.

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©2026 DJO, LLC
DJO FRANCE SAS
NOSS GmbH
DJO Singapore Pte Ltd
DJ ORTHO, LLC