OPTIMA CARE
MOLLITER

Fase Acuta / Acute Phase
Akutphase / Phase Aaigüe

MANUALE D'USO

INSTRUCTIONS FOR USE

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MODE D'EMPLOI

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DESCRIPTION OF THE DEVICES

The **Optima** line is indicated for the treatment of alterations in the morphology of the foot.

INTRODUCTION

Optima Molliter S.r.l. declares that Optima Molliter walkers are medical devices conceived and designed following the current state of the art, have been tested to be in compliance with the Regulation EU 93 2017/745, in its applicable sections. Optima Molliter walkers are in compliance with the International Guidelines of the IWGDF for the management of the diabetic foot ulcerations. https://iwgdfguidelines.org/guidelines/guidelines. Nevertheless, bear in mind that inappropriate use of product or underestimation of resulting dangers can cause injuries or harm. They must be correctly used to be effective. It is important to read this User Manual carefully before using the device and to report any serious incident occurring in relation to the device to the manufacturer and to the competent authority of the State in which the user and/or patient is located.

INTENDED USE

OPTIMA DIAB

- Diabetic foot ulcer.
- Foot ulcer.
- Rehabilitation after minor amputation (except for Lisfranc and Chopart amputations).
- Rehabilitation after tarsus surgery.
- Management of Charcot foot during transition from acute to chronic phase.
- Arthrosis.
- Malleolar fracture.
- Metatarsal fracture.

OPTIMA CLHEEL

- Management of heel ulcer (plantar and dorsal).
- Ulcer/laceration of Achilles tendon.
- Heel decubitus ulcer (plantar and dorsal).
- Postoperative rehabilitation of heel/Achilles tendon.

OPTIMA POST OP

- Diabetic foot ulcer with vascular complications in the leg, neuroischemic ulcer before and after revascularization (with additional use of the OPTIMA PUZZLE KIT 3X3).
- Management of the contralateral foot or management of patients at risk for foot ulceration.
- Metatarsal compound fracture.
- Rehabilitation after toe amputation.
- Preoperative stabilisation.

OPTIMA EUROPA

- Management of interdigital ulcer.
- Management after hallux valgus surgery.
- Management after hammer toe surgery.
- Post-operative management of the toes and Morton's neuroma.
- Fracture of toes.
- Nail excision.
- Plantar wart excision.
- Interdigital lesion.
- Toe lesions.

OPTIMA FREE

- Management of dorsal foot ulcer.
- After forefoot surgery.
- After toe surgery.
- After hallux valgus surgery.
- Recommended for containing functional bandages.

OPTIMA GROUND

- Metatarsal compound fracture.
- Rehabilitation after toe amputations.
- Preoperative stabilisation.
- Management of toes post-surgery.
- Toes fractures.

OPTIMA LIGHT 0.1

- Toes/feet fractures.
- Hallux Valgus.
- Hallux Rigidus.
- Hammer toes.
- Morton's Neuroma.
- Nail excision.
- Plantar wart excision.
- Interdigital lesion.
- Toe lesion.
- Post-surgery.
- Forefoot, midfoot and toe ulcers.

OPTIMA PUZZLE KIT 3x3 / OPTIMA KITCLHEEL

- Redistribuiting the plantar pression in case of pression ulcer of the forefoot, midfoot and heel and for the customized offloading of specific areas.
- It allows 27 different combinations of offloading.

OPTIMA PLTM

It is an insole that must be used inside the Optima Diab, Optima Post Op and Optima Ground braces, to protect the amputated foot and improve its position within the brace.

SIDE EFFECTS AND WARNINGS

No side effects related to the use of the devices were recorded.
It is advisable not to use the medical device directly on the skin, always cover the foot with bandages or sock. Skin lesions are possible due to excessive pressure of the skin in direct contact with the various parts of the device.
The product should be released and used only on medical prescription or on the advice of a licensed healthcare worker.
Do not use these devices if there is a decreased sensitivity resulting from anaesthetic administration. The 3x3 KIT multilayer insoles (standard on Optima Diab and Optima Clheel), the monolayer insole (standard on Optima Post Op, Optima Europa , Optima Free, Optima Ground and Optima Light) and the PLTM insole need to be handled exclusively by specialized operators (doctor or licensed healthcare worker).

Always check the integrity of the devices and notify the manufacturer in case of negative feedback.

- It is highly recommended to wear a contralateral protective shoe on the foot without lesions with the same height of the bottom sole of the device.
- Pay utmost care when walking if there are liquids or slippery substances on the floor.
- Do not expose the product to open flames or excess heat sources.
- Make sure that all fasteners are firmly secured before using the device to walk.

In the case of subjects with precarious walking, use appropriate medical aids to help walking such as crutches, tripods, walkers as indicated by your healthcare professional. Notify the Manufacturer in case of any adverse reaction which has not been mentioned in this manual. In the event of a serious accident, notify the Manufacturer and the Competent Authorities. To be used by subjects weighing up to 130 Kg.

IMPROPER USE

The manufacturer **Optima Molliter S.r.l.** shall not be held liable if the product is used under conditions listed below as improper use.
Failure to comply with the therapeutic indications / failure to comply with the rules and instructions of the manufacturer /non-compliance with or use contrary to current safety standards / significant alteration of dimensions and shape / significant alteration of the structure by changing-adding items / use product to practice sports, even if not extensive/ use of unsuitable cleaning products /driving automobiles while wearing products / unsuitable storage / device components are assembled/disassembled by non-qualified personnel / failure to comply with assembly or disassembly instructions / use the product like a tool or like a toy.

DURATION OF USE

The product retains its characteristics when properly used and stored over its life cycle (on average one season / 6 months), according to the environmental characteristics and its use. The OPTIMA KIT 3x3, OPTIMA KITCLHEEL and OPTIMA PLTM insoles are subject to greater wear and therefore their operational life is estimated to be around 21 days.

INSTRUCTIONS FOR USE

The product can be used directly or with the appropriate modifications to the insoles as expressly indicated. It is important not to modify the structure and the foreseen support points in order not to alter the functionality of the product.

CUSTOMIZATION OF THE OPTIMA PUZZLE KIT 3x3


The modular KIT 3x3 system included in Optima Diab and Optima Clheel and optional in Optima Post Op, Optima Europa, Optima Free, Optima Ground and Optima Light has been designed to redistribute plantar pressures and manage pressure ulcers / injuries of the forefoot, midfoot and heel.

WARNING!!!
The 3x3 KIT must be customized only by doctors, nurses, orthopedic technicians or qualified health personnel. **WARNING!!!** The 3x3 KIT must only be used within Optima Diab and Optima Clheel (1 kit is already series), Optima Post Op, Optima Europa, Optima Free, Optima Ground and Optima Light. (optional). Never use KIT 3x3 in normal footwear or other braces. Optima Molliter S.r.l. is not responsible for the use of KIT 3x3 in devices different from Optima Diab, Optima Clheel, Optima Post Op, Optima Europa, Optima Free, Optima Ground and Optima Light **WARNING!!!** The 3x3 KIT must always be used as a whole with all its components well connected to each other inside the brace.

KIT 3X3 consists of 9 interconnected parts in a Puzzle system forming 3 different small plantar insoles; the three plantar small insoles have 3 densities (**RED** = soft, **BEIGE** = medium, **BLUE**= rigid) and they can be combined vertically and horizontally. Plantar pressure can be easily redistributed by combining these insoles. The procedure consists in simple steps; to reduce pressure on an ulcerated/injured area, the KIT 3X3 can be customised by marking and cutting part of the plantar insole.

PATIENTS WEIGHING LESS THAN 90 KG

Mark the ulcerated/injured plantar point on the **BEIGE** insole
1.Cover the plantar injury/ulcer with a transparent sterilised film and using a lipstick or non-permanent marker mark the perimeter of the injury/ulcer directly over the transparent film (A). **2.** Insert the entire insole KIT 3X3 into the device and make sure that the **BEIGE** insole is in contact with the foot (**BLUE**-**BEIGE**-**RED**).
3. Put the patient's foot into the device and close the straps. Help the patient to walk (10 steps). In this way, the plantar surface will be marked at the injury point.
4. Now, open the device and remove the entire KIT 3X3.
5. Using a pen mark the perimeter of the imprint with a plus margin of about 6 mm around and 8 mm in front (B).
6. Using a scalpel with a 10 mm-blade perforate the perimeter of the mark with the blade at about 45° (C). **7.** Assemble the components according to the desired pressure redistribution (Diagram A) and insert the entire KIT 3X3 into the device.
8. Make sure that the insole with the hole is in the intermediate segment.
9. To avoid infections, always protect the ulcer/injury with suitable dressing before using the KIT 3X3. **10.** Once the Protective Sock has been put on, the device can now be used by the patient.


DIAGRAM A

PATIENTS WEIGHING MORE THAN 90 KG


Mark the ulcerated/injured plantar point on the **BLUE** insole.
1. Cover the plantar injury/ulcer with a transparent sterilised film and using a lipstick or non-permanent marker mark the perimeter of the injury/ulcer directly over the transparent film (D). **2.** Insert the entire insole KIT 3X3 into the device and make sure that the **BLUE** insole is in contact with the foot (**BLUE**-**BEIGE**-**RED**).
3. Put the patient's foot into the device and close the straps. Help the patient to walk (10 steps). In this way, the plantar surface will be marked at the injury point.
4. Now, open the device and remove the entire KIT 3X3.
5. Using a pen mark the perimeter of the imprint with a plus margin of about (E) 6 mm around and 8 mm in front. 6. Using a scalpel with a 10 mm-blade perforate on the **BLUE** insole the perimeter of the mark with the blade at about 45° (F).
7. Place the **BLUE** drilled part on the corresponding **BEIGE** part.
Mark the base of the conical hole with a pen. Drill a conical hole also on the **BEIGE** insole (G). **8.** Assemble the components according to the desired pressure redistribution (see Diagram B) and insert the entire KIT 3X3 into the device.
9. Make sure that the **BLUE** insole with the hole is in the intermediate segment and the **BEIGE** insole in the segment below.
10. To avoid infections, always protect the ulcer/injury with suitable dressing before using the KIT 3X3.
11. Once the Protective Sock has been put on, the device can be used by the patient.

DIAGRAM B

WARNING!!!
THE AREAS WITH THE OFF-LOADING HOLE MUST NEVER BE IN CONTACT WITH THE PLANTAR ULCER/INJURY. DO NOT USE THE DEVICE WITHOUT THE PLANTAR INSOLE. REPLACE THE ENTIRE KIT 3X3 EVERY 21 DAYS ON PRESCRIPTION OF A SPECIALIZED DOCTOR (IN THE CASE OF THE ULCER/LESION IS NOT HEALED).
For further information, request the video with the instructions for use to the following e-mail address: **info@optimamolliter.com** or scan the QR code.

PREPARATION OF THE PLTM PLANTAR INSOLE

1. Remove the adhesive film of the black Velcro band and attach it on the insole internal surface **2.** Attach the orange material with the Velcro section on the Velcro band **3.** Position the foot with minor amputation on the insole being careful that there is at least 1,5 cm between the foot and the orange material. **4.** Always protect foot and leg with adequate dressing. It is recommended to use a protective sock.

WARNING!!!
PLTM must be prescribed only by a specialist physician and customised only by physicians, carers, orthopaedic technicians or qualified healthcare personnel. PLTM must be used only inside Optima Diab, Optima Post Op and Optima Ground.
PLTM is designed for postoperative following toe amputation or transmetatarsal amputation. Do not use PLTM plantar insole in normal shoes or other devices. Optima Molliter is not responsable for the use of PLTM in different devices.

CLEANING OF THE PRODUCT

To clean the product, use warm or cold water and neutral soap; rinse once the operation is completed. Avoid any corrosive liquids such as alcohol, laundry detergents, bleach, ammonia, and solvents. Avoid using abrasive materials such as brushes, steel wool, and sand paper, or a knife blade to remove stains.Optima Molliter S.r.l. shall NOT be held liable for damages caused by abrasive products, by products which attack the materials, by direct exposure to heat sources.

DISPOSAL OF THE PRODUCT

Dispose of product responsibly. Refer to the local regulations on the disposal of inert products. Packaging can be disposed of in paper.

WARRANTY

The warranty refers only to materials used and to manufacturing defects and protects the device in accordance with the law only if recognized as defective by the Company. For the article, which is assumed to be defective, a complaint can only be requested by fulfill the "Complaints Form" which can be requested to the Company and sent to resi@optimamolliter.com only in the presence of:
- intact, unused article kept in its original packaging and joined by proof of purchase which must not be more than 10 days prior to the claim. With the prior authorization of Optima Molliter S.r.l., the item must be sent in the original packaging together with the proof of purchase, exclusively carriage paid and indicating the return authorization number communicated by Optima Molliter on the package. Returns that should arrive without affixing the authorization number clearly visible on the package will be rejected. Items used or shipped without compliance with the above specifications, or valued not covered by warranty, may be returned to the sender only at his charge. In the event of any recognized liability of the Company, the product will be returned by the Company to the reseller or customer.

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Optima Molliter Srl