#554

PAVIS S.p.A.

Sede / Headquarter:
Via G. Rossini, 11/B - 21020 Buguggiate (VA) - Italy
Tel. (+39) 0332 45 54 63
Fax (+39) 0332 45 97 48

www.pavis.it - info@pavis.it
pavispec@certimprese.it
R.I. / C.F. / P.IVA 01286640121
R.N. / F.C. / VAT IT01286640121

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
(UK)

1 - INDICATIONS
- Painful inflammatory/degenerative pathologies: Osteoporosis, osteomalacia, moderate osteolysis; back pain, curved back
- Post traumatic: After trauma of the dorso-lumbar-sacral spine with uniform tightness of the pelvis
- Post surgery: In stabilised outcomes of vertebroplasty, cytoplasty, discectomy

2 - CONTRA-INDICATIONS AND SIDE-EFFECTS
The choice between different models must be done according to the pathology and the necessities stated by the doctor. Do not wear it close to open flames. To avoid skin redness (caused by movements mechanical frictions) it is preferable to wear it over a cotton t-shirt.
Since this device is NON-INVASIVE and made with non-allergenic fabrics, no negative or side effects have ever been noted.
Any serious incident occurred in relation to the device use, should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

3 - MAINTENANCE AND CARE

This device is non-biodegradable. Even though it is made of natural fabrics, after use, avoid release the product to the environment. AT THE END OF ITS USABLE LIFE HAS TO BE DISCHARGED ALONG WITH SOLID ORDINARY WASTE.

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